

PRESS RELEASE

## “A clear view on flow”

Bronkhorst High-Tech proudly presents the BRIGHT Series, multifunctional Readout and Control units. These new R/C-modules offer the users of Bronkhorst Mass Flow Meters/Controllers the possibilities of local indication and operation, flowrate in actual values or percentage, fluid selection, setpoint, counter and alarms.

Further to these functionalities, the compact and stylish design distinguishes this comprehensive device from other products on the market. The BRIGHT models are available for lab-style instruments (see depicted image of an EL-FLOW<sup>®</sup> Mass Flow Controller) and for instruments with IP65 housing. Communication between such an instrument and the BRIGHT module is established over the RS-232 digital communication port.



The modern OLED display technology results in a bright and wide angle visibility. The BRIGHT module can be mounted either directly on the instrument, on an adjacent wall / panel or on the instrument's pipework. For the latter options we can offer an additional mounting bracket.

For more information contact:

**BRONKHORST HIGH-TECH B.V.**

Tel.: +31 573 458800

Fax: +31 573 458808

E-mail: [sales@bronkhorst.com](mailto:sales@bronkhorst.com)