


A066-GP99 - MESURE DE DÉBIT EN DÉPRESSION

APPLICATION NOTE A066-GP99

SUB-ATMOSPHERIC

Bronkhorst offers digital MASS-STREAM D-6300 mass flow meters and integrated flow controllers for gases with flow rates from 0,2 up to 5000 l/min Air equivalent.

The CTA (Constant Temperature Anemometry) inline measurement principle incorporates effectively an open pipe allowing mass flow measurement and control also with pressure ratings below 1 bar(a).



Application requirements


- $P_1 > \text{approximately } 300 \text{ mbar (a)}$
- Integrated valve control with optional valve function 'fails in place' in case of a power failure
- Gas flow suction through the instrument can be generated with a pump
- Process evaluation with a trial unit before final order placement

Industrial sectors

- Process
- Low pressure
- Vacuum
- Testing
- Research

Optional features

- Multi-functional displays
- Additional interface
- Up to 8 calibration curves in one instrument
- Degreased instruments for oxygen applications



Flow scheme

Recommended Products



MASS-STREAM D-6370 MFM

Débit min. 2...100 l/min
Débit max. 10...1000 l/min
Pression jusqu'à 20 bar
Boîtier robuste (IP65)
Option afficheur TFT intégré



MASS-STREAM D-6380 MFM

Débit min. 10...500 l/min
Débit max. 50...5000 l/min
Pression jusqu'à 20 bar
Boîtier robuste (IP65)
Option afficheur TFT intégré



BRONKHORST FRANCE S.A.S.

53 Rue Jacques Verniol

F-95370 Montigny-Les-Cormeilles (FR)

Tel. [+33 1 34 50 87 00](tel:+33134508700)

sales@bronkhorst.fr